Canine Distemper Virus Test Kit

* Principles

SensPERT Canine Distemper Virus Test Kit is designed to detect the antigens of canine distemper viruses in canine ocular discharge, nasal discharge. Two monoclonal antibodies in the kit specifically bind with different epitopes of the antigens. After being absorbed into the cellulose pad, the antigens of canine distemper virus move and bind with gold-colloid complex of monoclonal anti-canine distemper virus of the conjugate pad, forming Ab-Ag complex. This complex, then, forms Ab-Ag-Ab direct sandwich binding with the antibody of another monoclonal anti-canine distemper virus in the nitrocellulose membrane. The test results can be shown through the appearance of Control and Test lines where the principles of immunochromatography are used.

* Characteristics

- 1) One-step rapid test of canine distemper virus antigens
- 2) Easy examination with various samples
- 3) Rapid test results within 5 ~ 10 minutes
- 4) Expensive equipment not required
- 5) Easy storage and maintenance
- 6) High-purity and high-quality materials of the test kit increases its sensitivity and specificity.

* Materials (10 tests/kit)

- 1) Test 10 units
- 2) Diluent (buffer) 0.4 ml × 10 units
- 3) Disposable dropper 10 units
- 4) Swab for sample collection 10 units

* Composition

Specimen well (S : for dropping), Test line (T), and Control line (C) are marked on the device. Inside of it, the strip is composed of sample pad, conjugate pad, nitrocellulose membrane (test paper) and absorbent pad.

* Effect

Detection of canine distemper virus antigens from canine ocular discharge or nasal discharge.

* Uses

- 1) Specimen
Canine ocular discharge or nasal discharge
- 2) Collection and storage of specimen
 - a) 4 drops (100 μl) of with diluent is required.
 - b) Specimen can be stored at 2~8°C for 7 days.
 - c) Specimen should be kept at 22~25°C before use.
- 3) Test procedure
 - a) Use swab to collect canine tears (or eye mucus), or snivel and put the specimen into a sample tube containing 0.4 ml of diluent.
 - b) Swirl the diluent with the swab
 - c) Drop 4 drops (100 μl) into the specimen well (S)
 - d) Read the test results within 5 ~ 10 minutes

◆ Invalid result after 10 min

4) Interpretation of the results

Purple band should appear on the control line regardless of the test result. The presence of another band on the test line determines the result.

Control Line (C): The line should always appear regardless of the presence of the antigens of canine distemper virus. If this line does not appear, the test should

be considered invalid. This might be because of deterioration of the kit or the lack of specimen. It should be tested again with another kit.

Test Line (T): The presence of antigens of canine distemper virus determines the test line.

- 1) Negative: Control line only appears.

- 2) Positive: Both Test and Control lines appear.

- 3) Retesting
 - a) Both Test and Control lines do not appear.

- b) Test line only appears.

5) Further examinations

This test is for primary screening only. Consult veterinarians for further necessary examinations to obtain clinical test results.

* Precautions

- 1) Use for canine *in-vitro* diagnostic purposes only.
- 2) Use within 10 minutes after opening the pouch because the test kit is very sensitive to moisture and its effect may diminish.
- 3) Do not touch the result window.
- 4) Every specimen should be used with different droppers.
- 5) For testing, the buffer included should be used.
- 6) Do not use specimen being contaminated with microbes. It may turn out false positive or false negative.
- 7) Manage specimen carefully. It can deliver unknown virus or infectious bacteria.
- 8) Use disposable gloves when the infection of specimen is a possibility. And wash your hands later.
- 9) Dispose solid wastes after sterilizing them at 121°C for over 1 hour.
- 10) Do not use the kit when the pouch is torn, sealing is not good or expiration date is passed.

* Storage method and expiration date

The test kits stored at 2~30°C can be used for 24 months after manufacturing. Do not keep them in a refrigerator. However, if they are stored under cold circumstances, keep them at room temperature for 15 ~ 30 minutes before use.

* Exchange

The test kits are manufactured under strict quality control system. Nevertheless, if they are deteriorated during delivery, ask our distributors for exchange.

* Liability

The entire risk due to the performance of this product is assumed by the purchaser. The manufacturer shall not be liable for indirect, special or consequential damages of any kind resulting from use of this product.

Date issued : Sep. 08, 2016
A102-01E-03

VETALL
Manufactured by
VetAll Laboratories
#B-305, Samsung Techno Valley, 140, Tongil-ro,
Deogyang-gu, Gyeonggi-do, 10594, Korea
Tel : 82-2-22193459, Fax : 82-2-22193457

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