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**\* Introduction**

Canine parvovirus (CPV) and canine coronavirus (CCV) are potential pathogens for enteritis. Although their symptoms are quite the same, their virulence is different. CCV is the second leading viral cause of diarrhea in puppies with canine parvovirus being the leader. Unlike CPV, CCV infections are not generally associated with high death rates. CCV is not new to the canine population. Dual CCV-CPV infections were identified in 15-25% of cases of severe enteritis in the USA. Another study showed that CCV was found in 44% of fatal gastroenteritis cases that were initially identified as only CPV disease. CCV has been widespread among the canine population for many years. The age of dog is also important. If a disease occurs in puppy, it often leads to death. In mature dog the symptoms are more gentle. The possibility for healing is higher. Puppies less than twelve weeks of age are at the greatest risk and some especially weaker ones will die if exposed and infected. A combined infection leads to a much more severe disease than occurs with either CCV or CPV alone and is often fatal.

Group	Severity of Signs	Mortality Rate	Recovery Rate
CCV	+	0%	100%
CPV	+++	0%	100%
CCV + CPV	+++++	89%	11%

**\* Principles**

SensPERT™ CCV Ag / CPV Ag test kit is for the simultaneous detection of antigens of CCV and CPV in canine feces. The test results can appear on Control (C) and Test (T) lines where the principles of immunochromatography are used.

**\* Characteristics**

- 1) One-step rapid test of CCV Ag / CPV Ag
- 2) Rapid test results between 5 ~ 10 minutes
- 3) Expensive equipment not required
- 4) Easy storage and maintenance
- 5) High-purity and high-quality materials of the test kit increase its sensitivity and specificity.

**\* Materials (10 tests/kit)**

- 1) Test ..... 10 units
- 2) Diluent (buffer) ..... 1 ml x 10 unit
- 3) Disposable dropper for sample collection ..... 10 units
- 4) Swab for sample collection ..... 10 units

**\* Composition**

Specimen well (S : for dropping), Test line (T), and Control line (C) are marked on the device. Inside of it, the strip is composed of sample pad, conjugate pad, nitrocellulose membrane (test paper) and absorbent pad

**\* Effect**

Simultaneous detection of antigens of CCV and CPV in canine feces

**\* Uses**

- 1) Specimen  
Canine feces
- 2) Test procedure
  - a) Collect canine feces with a swab and put the swab into a sample tube containing 1 ml of buffer.  
**Do not use sticky or a lump of feces.**
  - b) Swirl the buffer with the swab.
  - c) Drop 4 drops (100 µl) into the specimen well (S).
  - d) Read the test results within 5~10 minutes.

♣ **Invalid result after 10 min**

**3) Interpretation of the results**

A purple band should appear on the control line regardless of the test result. The presence of another band on the test line determines the result.  
 Control Line (C): The line should always appear. If this line does not appear, the test should be considered invalid. This

might be because of impure buffer or the lack of specimen. It should be tested again with new buffer.  
 Test Line (T): The presence of CCV Ag and/or CPV Ag determines the test line.

1) CCV / CPV Negative

2) CCV Positive / CPV Negative

3) CCV Negative / CPV Positive

4) CCV / CPV Positive

5) Retesting

a) Both Test and Control lines do not appear.

b) Test line only appears.

**\* Further examinations**

This test is for primary screening only. Consult veterinarians for further necessary examinations to obtain clinical test results.

**\* Precautions**

- 1) Use for canine *in-vitro* diagnostic purposes only.
- 2) Use within 10 minutes after opening the pouch because the test kit is very sensitive to moisture and its effect may diminish.
- 3) Be careful of not touching the result window.
- 4) Every specimen should be used with different droppers.
- 5) For testing, the buffer included should be used.
- 6) Do not use specimen being contaminated with microbes, which may cause false positive or false negative result.
- 7) Deal with specimen carefully. They can deliver unknown virus or infectious bacteria.
- 8) Use disposable gloves when you suspect the infection caused by specimen. And wash your hands later.
- 9) Dispose solid wastes after sterilizing them at 121°C for over 1 hour.
- 10) Do not use the kit when its pouch is torn, sealing is not good in shape or expiration date is passed.

**\* Storage method**

The test kits stored at 2~30°C can be used for 24 months after manufacturing. Do not keep them in a refrigerator. However, if they are stored under cold circumstances, keep them at room temperature for 15 ~ 30 minutes before use.

**\* Exchange**

The test kits are manufactured under strict quality control system. Nevertheless, if they are deteriorated during delivery, ask our distributors for exchange.

**\* Liability**

The entire risk due to the performance of this product is assumed by the purchaser. The manufacturer shall not be liable for indirect, special or consequential damages of any kind resulting from use of this product.

Date issued : Sep. 08, 2016  
 A108-01E-02

**VETALL** Manufactured by  
**VetAll Laboratories**  
 #B-305, Samsong Techno Valley, 140, Tongil-ro,  
 Deogyang-gu, Goyang-si, Gyeonggi-do, 10594, Korea  
 Tel : 82-2-22193459, Fax : 82-2-22193457


# Instruction for CCV Ag / CPV Ag test


The amount of sample : 100 µl

## 1 Sampling

1. Do not use sticky or a lump of feces.
2. Do not collect too much feces.  
 Maximum quantity : 30 mg  
 (Fig : Swab with 30 mg of feces)


## 2 Use the Supernatant layer


## 3 Loading

4 drops (100 µl) of the diluted sample into each specimen well (S)  
**\*\*Do not overload the sample!!**

## 4 Reading

00:05:00~00:10:00

## 5 Discard

00:10:00

: Discard the used device  
**Invalid result after 10 min**