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



SensPERT[®] Giardia Ag Test Kit

Principles

SensPERT Giardia Test Kit is designed to detect the antigens of Giardia lamblia in feces. Two antibodies in the kit specifically bind with different epitopes of the antigens. After being absorbed into the cellulose pad, the antigens of Giardia lamblia move and bind with gold-colloid complex of anti-Giardia of the conjugate pad, forming Ab-Ag complex. This complex, then, forms Ab-Ag-Ab direct sandwich binding with the antibody of another monoclonal anti-Giardia in the nitrocellulose membrane. The test results can appear on Control and Test lines where the principles of immunochromatography are applied.

Characteristics

- 1) One-step rapid test of Giardia lamblia antigens
- 2) Rapid test results between 5 ~ 10 minutes
- 3) Expensive equipment not required
- 4) Easy storage and maintenance
- 5) High-purity and high-quality materials of the test kit increases its sensitivity and specificity.

Materials (10 tests/kit)

- 1) Test10 units
- 2) Diluent (buffer) 1ml x 10 units
- 3) Disposable dropper 10 units
- 4) Swabs for sample collection 10 units

Composition

Specimen well (S : for dropping), Test line (T), and Control line (C) are marked on the device. Inside it, the strip is composed of sample pad, conjugate pad, nitrocellulose membrane (test paper), and absorbent pad

Effect

Detection of Giardia lamblia antigens from canine or feline feces

Uses

- 1) Specimen
Canine or Feline feces
- 2) Storage of specimen
 - a) Store specimen at 2 ~ 8°C for use within 24 hours or at -25°C for later use.
 - b) Keep the specimen at 22 ~ 25°C before use.
- 3) Test procedure
 - a) Collect feces with a swab and put the swab into a sample tube containing 1ml of buffer.
Do not use sticky or a lump of feces.
 - b) Swirl the buffer with the swab.
 - c) Drop 4 drops (100µl) into the specimen well (S).
 - d) Read the test results within 5~10 minutes.

◆ Invalid after 10 minutes

- 4) Interpretation of the results
A purple band should appear on the control line regardless of the test result. The presence of another band on the test line determines the result.

Control Line (C): The line should always appear regardless of the presence of Giardia antigens. If this line does not appear, the test should be considered invalid. This is probably because of impure buffer and the lack of specimen. It should be tested again with new buffer.

Test Line (T): The presence of Giardia antigens determines the test line.

- 1) Negative: Control line only appears.

- 2) Positive: Both Test and Control lines appear.

- 3) Retesting
 - a) Both Test and Control lines do not appear.

- b) Test line only appears.

5) Further examinations

This test is for primary screening. Consult veterinarians for further necessary examinations to obtain clinical test results.

Precautions

- 1) Use for *In-vitro* diagnostic purposes only.
- 2) Use within 10 minutes after opening the pouch because the test kit is very sensitive to moisture and its effect may diminish.
- 3) Be careful of not touching the result window.
- 4) Every specimen should be used with different droppers.
- 5) For testing, buffer included only should be used.
- 6) Deal with specimen carefully. They can deliver unknown virus or infectious bacteria.
- 7) Use disposable gloves when you suspect the infection cause by specimen. And wash your hands later.
- 8) Dispose solid wastes after sterilizing them at 121°C for over 1 hour.
- 9) Do not use the kit when its pouch is torn, sealing is not good shape and expiration date passes.

Storage method

The test kits store at 2 ~ 30°C can be used for 24 months after manufacturing. Do not keep them in a refrigerator. However, if they are stored under cold circumstances, keep them at room temperature for 15 ~ 30 minutes before use.

Exchange

The test kits are manufactured under strict quality control system. Nevertheless, if they deteriorate during delivery, ask our distributors for exchange.

Liability

The entire risk due to the performance of this product is assumed by the purchaser. The manufacturer shall not be liable for indirect, special or consequential damages of any kind resulting from use of this product.

Date issued : Feb 16, 2017
A113-01E-01

Manufactured by
VetAll Laboratories
#B-305 Samsung Techno Valley, 140 Tongil-ro, Deogyang-gu,
Goyang-si, Gyeonggi-do, 10594, Korea
TEL. +82-2-2219-3459 FAX. +82-2-2219-3457


Instruction for Giardia Ag test

The amount of sample : 100 µl

1 Sampling

1. Do not use sticky or a lump of feces.
2. Don't collect too much feces.
Maximum quantity : 30 mg
(Fig : Swab with 30 mg of feces)

Swirl the buffer with the swab

2 Use the Supernatant layer


Wait until feces particles sink to the bottom


Supernatant layer for use

Feces particle layer


3 Loading

4Drops (100 µl) of the diluted sample
****Don't overload the sample!!**

4 Reading

00:05:00~00:10:00

Read the result between 5~10 min

5 Discard

00:10:00

: Discard the used device
Invalid result after 10min